

WDS

Digital X-ray sensor
Bluetooth® Technology


Minimum System Requirements	
Operating System	Microsoft® Windows® VISTA or XP Service Pack 2/3 Apple® Mac OS X 10.5 Leopard
Ports	Integrated Bluetooth® or USB 2.0
Processor	1.0 GHz or faster
Memory	512 MB
Display Setting	1024x768 at 32 bit true color
Technical Specifications	
Sensor:	CCD + FOP + CsI
Pixel Size:	20x20µm
Intrinsic resolution:	CCD 25 lp/mm
Sensitivity:	30mR half well
Dynamic range:	12 bit, 4096 gray shades
Size 1	
Pixel Matrix:	1000 x 1500
Active Area:	20 x 30 mm
External Dimensions:	26.4 x 38.8 x 5.3 mm
Size 2	
Pixel Matrix:	1200 x 1700
Active Area:	24 x 34 mm
External Dimensions:	30.7 x 42.7 x 5.3 mm

MWDSGB91S00

2009/10/15

Data subject to change without notice





Three-layer Sensor Technology

State-of-the-art X-ray image acquisition. Made of three different layers encapsulated in a protective shell, each layer contributes to final image quality.

CsI

This Caesium Iodide scintillator is the first one to intercept the X-ray beam converting it into visible light. It is manufactured with a vertical growth process that generates a columnar microstructure, which grants unsurpassed image quality.

FOP

This Fibre Optics Plate protects the sensor from direct X-ray penetration, allowing years of use without image deterioration. Moreover, the vertical fibres preserve the image resolution while the light propagates through the three layers in the sensor.

CCD

This is the acquisition device. It converts the light into a digital image. With perfectly square acquisition cells of 20µm, it can reach the maximum resolution of 25 lp/mm with 12bit encoding and 4096 gray shades, way above your intra-oral imaging needs. So you can be sure no detail is missed.

Built with your patient's health in mind

When it comes to diagnosis, images are not enough. It is essential that those images are accurate and clear.


WDS features a highly sophisticated X-ray sensor, capable of image acquisition at minimal exposure dosage, among the lowest in the industry. So you can control any step of dental treatment without increasing radiation. More control, better diagnosis, lower dosage. All in all, it's all about care.

WDS makes it all simple to achieve. As easy to handle as a cell phone and as simple to use as a remote control, this little device has a big technological heart. No wires, no strings, no power supply.



WDS is a pocket-size handheld device, powered by normal rechargeable batteries. Carry it from room to room freely, saving time and costs.



Sensor with FOP



Sensor without FOP



Bluetooth® Technology

Fast and reliable, this widespread industry standard has become very popular on mobile phones.

The small Bluetooth® receiving antenna plugs directly into the USB port of your PC.

No need for expensive control boxes and receivers, nor additional investments to share the system among multiple treatment rooms. Just a small Bluetooth® antenna plugged into each computer.



- Sensor name
- Transmission status
- Computer name
- Battery level

The Digital Way

Save time and costs, reduce X-ray exposure, improve your diagnostics.

Enhance, share, store your X-ray images easily, thanks to the treatment software. The WDS sensor works with existing X-ray generators. It is available in two sizes, so you can choose the one that better suits your needs. Specific aiming devices are designed to facilitate sensor positioning.

WDS is intraoral X-ray technology made simple, so nothing else stands in the way between you and the digital choice. Going digital is no longer an option. Better workflow. Better care.

